

Gene Name	Valosin Containing Protein (VCP)
Protein Name	
Synonyms	CDC48, TERA, p97
Uniprot ID	<u>P55072</u>
Description	The VCP gene encodes a member of the AAA ATPase family of proteins. The encoded protein plays a role in protein degradation, intracellular membrane fusion, DNA repair and replication, DNA damage, regulation of the cell cycle, autophagy, transport, activation of the NF-kappa B pathway and the Ubl conjugation pathway [1-3]. VCP has been implicated in human diseases such as: amyotrophic lateral sclerosis (ALS), Charcot-Marie-Tooth disease, inclusion body myopathy with paget disease of bone and frontotemporal dementia (IBMPFD), neurodegeneration, and neuropathy [4-8].
Antigen Coverage	Full length construct, 806 aa
Antigen Sequence (Antigen - <u>C-terminal</u> linker and Avi tag)	SMASGADSKGDDLSTAILKQKNRPNRLIVDEAINEDNSVVSLS QPKMDELQLFRGDTVLLKGKKRREAVCIVLSDDTCSDEKIRM NRVVRNNLRVRLGDVISIQPCPDVKYGKRIHVLPIDDTVEGITG NLFEVYLKPYFLEAYRPIRKGDIFLVRGGMRAVEFKVVETDPS PYCIVAPDTVIHCEGEPIKREDEEESLNEVGYDDIGGCRKQLA QIKEMVELPLRHPALFKAIGVKPPRGILLYGPPGTGKTLIARAV ANETGAFFFLINGPEIMSKLAGESESNLRKAFEEAEKNAPAIIFI DELDAIAPKREKTHGEVERRIVSQLLTLMDGLKQRAHVIVMAA TNRPNSIDPALRRFGRFDREVDIGIPDATGRLEILQIHTKNMKL ADDVDLEQVANETHGHVGADLAALCSEAALQAIRKKMDLIDLE DETIDAEVMNSLAVTMDDFRWALSQSNPSALRETVVEVPQVT WEDIGGLEDVKRELQELVQYPVEHPDKFLKFGMTPSKGVLFY GPPGCGKTLLAKAIANECQANFISIKGPELLTMWFGESEANVR EIFDKARQAAPCVLFFDELDSIAKARGGNIGDGGGAADRVINQ ILTEMDGMSTKKNVFIIGATNRPDIIDPAILRPGRLDQLIYIPLPD EKSRVAILKANLRKSPVAKDVDLEFLAKMTNGFSGADLTEICQ RACKLAIRESIESEIRRERERQTNPSAMEVEEDDPVPEIRRDH FEEAMRFARRSVSDNDIRKYEMFAQTLQQSRGFGSFRFPSG NQGGAGPSQGSGGGTGGSVYTEDNDDDLYG <u>SSKGGYGLND IFEAQKIEWHE</u>
Z-VCP-10 scFv Sequence (V _H - <u>linker</u> - V _L)	EVQLLESGGGLVQPGGSLRLSCAASGFTFSSYAMSWVRQAP GKGLEWVSAISGSGGSTYYADSVKGRFTISRDNSKNTLYLQM NSLRAEDTAVYYCARSYVGTMDYWGQGTLVTVSS <u>GGGGSG</u> <u>GGGSGGGGS</u> DIQMTQSPSSLSASVGDRVTITCRASQSISSYL NWYQQKPGKAPKLLIYAASSLQSGVPSRFSGSGSGTDFTLTI SSLQPEDFATYYCQQSPGYYLSTFGQGTKLEIK
Antibody Fragment Generation and Production	Z-VCP-10 was generated by phage display technology from a human synthetic scFv library. The scFv was produced by bacterial expression in <i>E. coli</i> , with secretion to the periplasm, and purified by protein A affinity. The final protein batch of Z-VCP-10 scFv also carries a C-terminal His ₆ tag and a 3xFLAG tag.
Validated Applications	ELISA, HTRF, Luminex, SPR, and IP-MS.



Enzyme-Linked ImmunoSorbent Assay (ELISA):

Method Description	ELISA screens were performed to assess binding of selected single chain fragment variable (scFv) antibody clones. Each clone was tested in duplicate against its target antigen and a streptavidin control. ELISA screens were performed in 384 well plates coated with [1ug/mL] streptavidin. Biotinylated antigens were diluted in PBT buffer (0.5% BSA, 0.05% Tween 20 in PBS) to [1ug/mL] and added prior to the addition of scFv bacterial culture supernatants. Bound scFvs, which contain a FLAG-tag, were detected using an HRP conjugated anti-FLAG M2 antibody (Sigma, [1:10,000]). TMB substrate solution (Thermo Scientific) was used as a chromogenic substrate, before the final addition of 1M H ₂ SO ₄ stop solution. Plates were washed four times in between each step listed above with PBS-T (0.05% Tween 20 in PBS). Absorbance was measured at 450nm. Only scFv clones with a high enough absorbance ratio (Target / Streptavidin) were selected for further analysis. ELISA results for 7/CP-10 are
	selected for further analysis. ELISA results for Z-VCP-10 are shown below (blank indicates no scFv was added).

Z-VCP-10 binding in ELISA screen:



Homogenous Time Resolved Fluorescence (HTRF):

Method Description	HTRF technology was employed to robustly and reliably
	determine in solution the binding of each scFy antibody to its
	determine, in solution, the binding of each set v antibody to its
	target antigen. Sequence unique clones, that tested positive in
	ELISA, were evaluated by HTRF. Supernatants from bacterial
	cultures containing scFvs were diluted in assay buffer (0.1% BSA
	in PBS), then mixed with biotinylated target antigens to a final
	concentration of 50nM. Donor anti-FLAG M2 antibody labeled
	with Terbium (Cisbio) and acceptor Streptavidin labeled with
	XL665 (Cisbio) were added to each scFv in 384 well plates.
	Samples were incubated for 2 hours at RT. The binding signal
	(665nm) and the background/noise signal (615nm) were
	measured using a EnVision instrument (PerkinElmer). HTRF
	ratios (665nm Acceptor / 620nm Donor) were calculated. Delta



ratios (ΔR) were determined by subtracting the background from
the mean emission fluorescence ratio for each sample. Delta
ratios were plotted to compare the energy transfer signal for each
scFv to its non-target control. Results for Z-VCP-10 are shown in
the graph below (blank indicates no scFv was added).



Luminex Assay:

Method Description	The specificity of the anti-VCP scFv antibodies were evaluated against 48 different antigens (including VCP) using a bead based multiplex Luminex assay. Color-coded magnetic beads (Magplex Luminex Corp.) were activated with [5mg/mL] NHS and EDC in activation buffer (0.1M NaH ₂ PO ₄ , pH 6.2), then added to 96 well plates and incubated with shake at RT for 20min. Beads were washed twice with MES buffer (0.05M 2-(N-morpholino) ethane sulfonic acid, pH 5.0) before adding neutravidin and incubating for 2 hours at RT with shake. Beads were then washed twice with PBS-T before storage buffer (Blocking Reagent for ELISA (Roche) with 0.1% ProClin300 (Sigma) was added and incubated overnight at 4°C. Storage buffer was removed and biotinylated proteins (1µM in PBS) were added to neutravidin coupled beads then incubated for 1 hour at RT with shake. Plates were washed 3x with PBS-T, storage buffer was added, and incubated at 4°C overnight. scFv antibodies were diluted in assay buffer (3% BSA, 0.05% Tween 20, and [10µL/mL] neutravidin in PBS) and incubated for 1 hour at 4°C. A beadstock was prepared with coupled beads in storage buffer at a concentration of 100beads/ID/µL. The beadstock solution (5µLs) was mixed with each scFv (45µL) in 384 well plates and incubated for 1 hour at RT with shake. Plates were analyzed on a FLEXMAP3D (Luminex Corp). Luminex results for Z-VCP-10 are plotted below. Graph illustrates the degree of specificity of Z-VCP-10 to bind to its respective antigen (VCP) and none of the
	other 47 analytics (blank indicates no SCFV was added).



Luminex results for Z-VCP-10:



Z-VCP-10 Blank

Surface Plasmon Resonance (SPR) Affinity Measurements:

Method Description	Kinetic measurements were performed using SPR on a Biacore T200 (GE Healthcare). The anti-FLAG M2 antibody (Sigma) was used to capture the FLAG-tagged scFv before antigen injection. Single Cycle Kinetics was performed. Accordingly, five concentrations of antigen are sequentially injected in increasing order within the same cycle before regeneration with Glycine-HCl, pH 2.5. Antigen concentrations tested against Z-VCP-10 include: 0.16, 0.8, 4, 20, and 100nM with a flowrate of 30µL/min. The
	kinetic constants were calculated using the Biacore T200 Evaluation Software 3.1 and the 1:1 Langmuir binding model, after removal of the reference cycle (running buffer instead of antigen) and the reference channel. Transient spikes from <i>e.g.</i> air bubbles were also removed from the obtained sensorgrams. The Single Cycle Kinetic sensorgram for Z-VCP-10 is shown below.

Single Cycle Kinetics curve for Z-VCP-10:



Susanne Gräslund, Structural Genomics Consortium Susanne.Gräslund@ki.se



Immunoprecipitation – Mass Spectrometry (IP-MS):

Method Description	Immunoprecipitation followed by mass spectrometry was
	employed to verify the interaction specificity of each scFv
	antibody with its intended target. Immunoprecipitation was
	performed on HEK293 cell lysates expressing endogenous
	VCP/TERA protein as previously described [9]. Obtained samples
	were loaded onto a Dionex Ultimate 3000 HPLC system (Thermo
	Fisher Scientific) coupled to an o-QTOF impact II™ (Bruker
	Daltonics) mass spectrometer with a Captive Spray ion source
	(Bruker Daltonics). Data was analyzed with Protein Scape
	software (Bruker Daltonics) using Mascot search engine (Matrix
	Science Ltd.). Normalized spectral abundance factor (NSAF)
	values were calculated as previously described [9].

Z-VCP-10 was able to capture endogenous VCP/TERA from HEK293 cell lysates. The target protein was at the top of the list of immunoprecipitated proteins with an NSAF value of 845 and 32 unique peptides detected. No other specific interactors were found. More detailed data can be shared upon request.

Comments and Contact:

Two additional scFv antibodies against VCP, from the same selection campaign, passed all validation criteria, but with slightly lower values. Sequences for those can be provided upon request. For all inquiries, please contact Dr. Susanne Gräslund (<u>susanne.graslund@ki.se</u>) at SGC Karolinska.

References:

- 1. Acs, K., et al., *The AAA-ATPase VCP/p97 promotes 53BP1 recruitment by removing L3MBTL1 from DNA double-strand breaks.* Nat Struct Mol Biol, 2011. **18**(12): p. 1345-50.
- 2. Cloutier, P., et al., A newly uncovered group of distantly related lysine methyltransferases preferentially interact with molecular chaperones to regulate their activity. PLoS Genet, 2013. 9(1): p. e1003210.
- 3. Meerang, M., et al., *The ubiquitin-selective segregase VCP/p97 orchestrates the response to DNA double-strand breaks.* Nat Cell Biol, 2011. **13**(11): p. 1376-82.
- 4. Gonzalez, M.A., et al., A novel mutation in VCP causes Charcot-Marie-Tooth Type 2 disease. Brain, 2014. 137(Pt 11): p. 2897-902.
- 5. Jerath, N.U., et al., *Rare Manifestation of a c.290 C>T, p.Gly97Glu VCP Mutation.* Case Rep Genet, 2015. **2015**: p. 239167.
- Johnson, J.O., et al., Exome sequencing reveals VCP mutations as a cause of familial ALS. Neuron, 2010. 68(5): p. 857-64.
- 7. Watts, G.D., et al., Novel VCP mutations in inclusion body myopathy associated with Paget disease of bone and frontotemporal dementia. Clin Genet, 2007. **72**(5): p. 420-6.
- 8. Watts, G.D., et al., Inclusion body myopathy associated with Paget disease of bone and frontotemporal dementia is caused by mutant valosin-containing protein. Nat Genet, 2004. **36**(4): p. 377-81.
- Persson, H., et al., Antibody Validation by Immunoprecipitation Followed by Mass Spectrometry Analysis, in Synthetic Antibodies: Methods and Protocols, T. Tiller, Editor. 2017, Springer New York: New York, NY. p. 175-187.